





► JY-DS1 & JY-LC170

Features:

- 1:Classic style
- 2:Lightweight and easy to use
- 3,Option:2.5X/3.5X

JY-DS2 & JY-LC180-B ◀

Features:

- 1:Classic loupes
- 2:Spotlight
- 3,Option:2.5X/3.5X



► JY-MGB & JY-LF180-B

Features:

- 1:Can Self-matching Private lenses
- 2:spotlight
- 3,Option:2.5X/3.5X

JY-ELD & JY-LC180-B ◀

Features:

- 1:Loghter distance(520mm)
- 2:Spotlight
- 3,Option:2.5X/3.5X





► JY-MFB & JY-LF180

Features:

- 1, Metal frame
- 2, Spotlight
- 3, Option: 2.5X/3.5X

JY-DS2 & JY-LF200 ◀

Features:

- 1, Wireless lights
- 2, Option: 2.5X/3.5X



► JY-FDG

Features:

- 1, High multiples
- 2, Option: 5X/6X
- 3, Working distance:
400, 450, 500, 550mm (Option)

JY-EKJ ◀

Features:

- 1, Compact enough to clip directly onto glasses
- 2, Option: 2.5X/3.5X





▶ **JY-LF200-A**
Wireless lights (Spot light)



▶ **JY-LF200-B**
Wireless lights



▶ **JY-LC170**
(Clip option)



▶ **JY-LF180-A** (Spot light)
(Clip option)



▶ **JY-LF180-B** (Spot light)
(Clip option)



▶ **JY-LF180-C** (Spot light)
(Adjust the size)
(Clip option)




► JY-ES2 & JY-LF180-A

Features:

- 1:Classic headband loupes
- 2:Spotlight
- 3,Option:2.5X/3.5X/5X/6X



► JY-AI02 & JY-LF190-A (Option:2.5X/3.5X)



► JY-AI02 & JY-LF190-B (Option:2.5X/3.5X)



► JY-ES2-B (Option:2.5X/3.5X)



► JY-ES2-C(Option:5X/6X) Working distance:400, 450, 500, 550mm(Optional)

► JY-PSL-A



Lens option: 2.5X / 3.5X



Magnifying glass



Filters



Charger (US/EU)




Charging base



Battery

► JY-PSL-B



Lens option: 5 X / 6X



Magnifying glass



Filters



Charger (US/EU)



Charging base



Battery

NEW

BLOCKBUSTER

NEW PRODUCT

NEW PRODUCT LAUNCH

NEW GENERATION HEADLIGHTS

NEW RESEARCH AND DEVELOPMENT/NEW PRODUCTS/NEW QUALITY



[NEW PRODUCTS ARE COMING]

Product Introduction

- [Integrated design] Convenient to carry.
- [Head mounted design] Convenient to use.
- [Adjustable headband] Suitable for more people.
- [Adjustable Brightness] Suitable for different usage scenarios.
- [Two colors] Available in black/white.



Light effect

Specifications

Product name:	LED bulb:
Integrated headlight	3W
Lifetime of led:	Battery:
10000h	1300mAh
Input voltage:	Color temperature:
AC100-240V/50-60Hz	5000k
Dominant wavelength:	0 cm distance luminosity:
487nm	480000lux(Max power)
10 cm distance luminosity:	20 cm distance luminosity:
95000lux	25000lux
30 cm distance luminosity:	40 cm distance luminosity:
10000lux	7000lux
Cotinues working time (highest intensity):	
≥2.5h	
Product net weight:	Packing size:
About 192g	25.4*22.4*12.8cm

NEW

BLOCKBUSTER

NEW PRODUCT

NEW PRODUCT LAUNCH

NEW SECOND-GENERATION HEADLIGHTS


NEW RESEARCH AND DEVELOPMENT/NEW PRODUCTS/NEW QUALITY



[NEW PRODUCTS ARE COMING]

Product Introduction

- [Integrated design] Convenient to carry.
- [Head mounted design] Convenient to use.
- [Adjustable headband] Suitable for more people.
- [Adjustable Brightness] Suitable for different usage scenarios.
- [Two colors] Available in black/white.



Light effect

Specifications

Product name:	LED bulb:
Integrated headlight	5W
Lifetime of led:	Battery:
10000h	1300mAh
Input voltage:	Color temperature:
AC100-240V/50-60Hz	5000k
Dominant wavelength:	0 cm distance luminosity:
487nm	300000lux(Max power)
10 cm spot size:	20 cm spot size:
4.5cm	8cm
30 cm spot size:	40 cm spot size:
10cm	14cm
Cotinues working time (highest intensity):	
≥2h	
Product net weight:	Packing size:
About 205g	25.4*22.4*12.8cm

NEW

BLOCKBUSTER

NEW PRODUCT

NEW PRODUCT LAUNCH

NEW THIRD-GENERATION HEADLIGHTS


NEW RESEARCH AND DEVELOPMENT/NEW PRODUCTS/NEW QUALITY



[NEW PRODUCTS ARE COMING]

Product Introduction

- [Integrated design] Convenient to carry.
- [Head mounted design] Convenient to use.
- [Adjustable headband] Suitable for more people.
- [Adjustable Spot] Meeting Different User Needs.
- [Adjustable Brightness] Suitable for different usage scenarios.
- [Two colors] Available in black/white.



Light effect

Specifications

Product name:	LED bulb:
Integrated headlight	5W
Lifetime of led:	Battery:
10000h	1300mAh
Input voltage:	Color temperature:
AC100-240V/50-60Hz	5000k
Dominant wavelength:	0 cm distance luminosity:
487nm	300000lux(Max power)
10 cm spot size:	20 cm spot size:
0.5-2.5cm	0.5-4cm
30 cm spot size:	40 cm spot size:
0.5-5cm	1-9cm
Cotinues working time (highest intensity):	
≥2h	
Product net weight:	Packing size:
About 205g	25.4*22.4*12.8cm

NEW

BLOCKBUSTER

NEW PRODUCT

NEW PRODUCT LAUNCH

EMBEDDED MAGNIFYING GLASS

NEW RESEARCH AND DEVELOPMENT/NEW PRODUCTS/NEW QUALITY



[NEW PRODUCTS ARE COMING]

Product Introduction

- [Embedded Design] Sensorless Wear.
- [Optical Glass Lens] The lens adopts German Schott optical lens, which provides clear imaging in the full field of view.
- [Optional Pupil Distance] Suitable for more people.

Optional Pupillary Distance

56 / 58 / 60 / 62 mm

Specifications

Magnification	Work distance	Field of view	Depth of vision
2X	320-420mm	100mm	90mm

NEW

BLOCKBUSTER

NEW PRODUCT

NEW PRODUCT LAUNCH

CORNER MAGNIFYING GLASS

NEW RESEARCH AND DEVELOPMENT/NEW PRODUCTS/NEW QUALITY



[NEW PRODUCTS ARE COMING]

Product Introduction

- [Glasses style design] Easy to use, seamless wearing.
- [Optical Glass Lens] The lens adopts German Schott optical lens, which provides clear imaging in the full field of view.
- [Corner Design] Unique optical corner structure design, bidding farewell to bending down and enjoying a comfortable working experience.

Specifications

Magnification	Work distance	Field of view	Peptth of visio
4X	400-500mm	50mm	50mm

NEW

BLOCKBUSTER

NEW PRODUCT

NEW PRODUCT LAUNCH

CORNER MAGNIFYING GLASS WITH LIGHT

NEW RESEARCH AND DEVELOPMENT/NEW PRODUCTS/NEW QUALITY



[NEW PRODUCTS ARE COMING]

Product Introduction

- [Corner Design] Unique Optical Corner Structure Design, Don't lower your head and bend down, enjoy a comfortable working experience.
- [Optical Glass Lens] The lens adopts Schott Optics from Germany Lens, clear full field imaging.
- [Paired with Soft Light Source] Paired with Lightweight High Brightness Spotlight Headlights Ultra long battery life.

Specifications(Corner Loupe)

Magnification	Work distance	Field of view	Peptth of visio
4X	400-500mm	50mm	50mm

Specifications(Headlight)

LED blud:	Battery:
5W	3500mAh
Lifetime of led:	Adpter voltage:
10000h	5V / 1A
Input voltage:	Color temperature:
AC100-240V/50-60Hz	5000k
Light intensity:	0 cm distance luminosity:
4 levels	200000lux(Max power)
10 cm spot size:	20 cm spot size:
3cm	6.5cm
30 cm spot size:	40 cm spot size:
9cm	12cm
Cotinues working time (highest intensity):	
≥7h	

NEW

BLOCKBUSTER

NEW PRODUCT

NEW PRODUCT LAUNCH

CORNER MAGNIFYING GLASS

NEW RESEARCH AND DEVELOPMENT/NEW PRODUCTS/NEW QUALITY



[NEW PRODUCTS ARE COMING]

Product Introduction

- [Glasses style design] Easy to use, seamless wearing.
- [Optical Glass Lens] The lens adopts German Schott optical lens, which provides clear imaging in the full field of view.
- [Corner Design] Unique optical corner structure design, bidding farewell to bending down and enjoying a comfortable working experience.

Specifications

Magnification	Work distance	Field of view	Peptth of visio
4X	400-500mm	50mm	50mm

NEW

BLOCKBUSTER

NEW PRODUCT

NEW PRODUCT LAUNCH

CORNER MAGNIFYING GLASS WITH LIGHT

NEW RESEARCH AND DEVELOPMENT/NEW PRODUCTS/NEW QUALITY



[NEW PRODUCTS ARE COMING]

Product Introduction

- [Corner Design] Unique Optical Corner Structure Design, Don't lower your head and bend down, enjoy a comfortable working experience.
- [Optical Glass Lens] The lens adopts Schott Optics from Germany Lens, clear full field imaging.
- [Paired with Soft Light Source] Paired with Lightweight High Brightness Spotlight Headlights Ultra long battery life.

Specifications(Corner Loupe)

Magnification	Work distance	Field of view	Pept of visio
4X	400-500mm	50mm	50mm

Specifications(Headlight)

LED blud:	Battery:
5W	3500mAh
Lifetime of led:	Adpter voltage:
10000h	5V / 1A
Input voltage:	Color temperature:
AC100-240V/50-60Hz	5000k
Light intensity:	0 cm distance luminosity:
4 levels	200000lux(Max power)
10 cm spot size:	20 cm spot size:
3cm	6.5cm
30 cm spot size:	40 cm spot size:
9cm	12cm
Cotinues working time (highest intensity):	
≥7h	